

Pathogenicity Calculator 3.00.xx

Installation instruction

Prerequisites

- Server with Genboree, version 2.06.00 or newer
- Access to console on the Genboree server as user “genboree”
- Access to the Allele Registry server, version 0.12.xx or later
- The following group will be created in Genboree:
 - pcalc_resources
 - pcalc_cache
- The following knowledge bases will be created in GenboreeKB:
 - pcalc_resources
 - pcalc_cache

Installation

1. Login to console on Genboree server as user “genboree”.
2. Download and decompress installation package with Pathogenicity Calculator.
3. Enter to the main package directory and run script *install.sh*.
4. Users in Pathogenicity Calculator may be created by using Genboree Redmine plugin “register_clingen_user”. It must be installed from Genboree Redmine and comes with Genboree.